

KINEMATRACE AI — CLINICAL GAIT ASSESSMENT REPORT

Markerless Kinematic Motion Analysis & Multi-Agent Risk Screening

1. PATIENT INFORMATION

Patient ID:	KT-2026-P902	Assessment Type:	Outpatient Gait Screening
Age Group:	2	Source Video:	gait_scan.mp4

2. AGENT 1 — VIDEO QUALITY VALIDATION

Validation Status:	PASS (Score: 100/100)
Landmark Tracking:	33/33 MediaPipe 3D Pose Keypoints Tracked
Camera Evaluation:	Satisfactory resolution, lighting, and frame rate

3. AGENT 2 — GAIT KINEMATIC MEASUREMENTS

Parameter	Patient Measurement	Normative Reference	Deviation Status
Gait Symmetry Index	83.96%	≥ 85.0%	Asymmetric
Mean Asymmetry Index	16.0%	≤ 15.0%	Elevated Asymmetry
Left Knee ROM	61.5°	110.0° - 140.0°	Extracted
Right Knee ROM	58.2°	110.0° - 140.0°	Extracted
Bilateral ROM Deficit	3.3°	≤ 5.0°	Normal
Hip Flexion ROM	120.0°	120.0° - 125.0°	On Benchmark

4. AGENT 3 — CLINICAL RISK ASSESSMENT

Screening Risk Level:	HIGH RISK (Severity: SIGNIFICANT)
Affected Limb:	RIGHT Limb
Explainable Reasoning:	Gait analysis measured mean asymmetry of 16.0% and bilateral ROM deficit of 3.3°.
Clinical Recommendation:	Schedule follow-up pediatric physical therapy assessment.

MEDICAL SAFETY DISCLAIMER: This automated gait screening report is generated by KinemaTrace AI for clinical decision support and markerless motion tracking. It is not a formal medical diagnosis. All findings and kinematic metrics must be interpreted by a qualified healthcare professional.